

Quantitative sulphur determination in BaSO₄ with rapid CS cube

BaSO₄ has a decomposition temperature of 1580°C and the remaining BaO can retain the release of SO₂. A quantitative sulphur analysis therefore requires a very high temperature during combustion and binding of the formed BaO.

Task

| | Instrument | Sample |
|-------------------|----------------------|----------------------------|
| Basic instrument: | rapid CS cube | Amount: 10 mg |
| Mode: | S | Consistency: powder |
| Peripheral: | manual sample former | Preparation: not necessary |

Specification

The sample is weighed in tin foil and formed to a pellet by means of the manual sample former after adding WO₃ powder 3 times of the sample weight to bind the BaO during combustion.

Procedure

| Sample no. BaSO ₄ | S [%] |
|------------------------------|--------------------------------|
| 1 | 13.80 |
| 2 | 13.69 |
| 3 | 13.59 |
| 4 | 13.86 |
| 5 | 13.71 |
| 6 | 13.63 |
| Mean | 13.71 ± 0.10 (theory 13.72% S) |

Results

The 100% S recovery and excellent precision show the ability to decompose compounds at an actual temperature of 1580°C and more due to tin oxidation during combustion.