

The sulphur content of Dried Distillers Grains measured with the rapid CS cube

Task

Dried distillers grains (DDG) are a by-product of the distillation process during spirit and bio-fuel production. Distillers grains are an excellent digestible protein and energy source for livestock, especially ruminants. For quality control the sulphur content of DDG is routinely determined.

Instrument	Sample
Basis: rapid CS cube	Quantity: 80-100 mg
Mode: S	Consistency: particulates 1-3 mm
Periphery: micro balance (readability 0.01 mg)	Preparation: not necessary

Specification

The samples are weighed into tin boats with an additive of WO₃. For calibration and the daily factor determination sulphanilamide is used.

Procedure

Each sample has been analysed four times. The average S content of six different dried distillers grains samples including their absolute standard deviations are presented below.

Sample	S [%]
DDG-1	1.064 ± 0.028
DDG-2	1.051 ± 0.037
DDG-3	0.989 ± 0.049
DDG-4	1.057 ± 0.040
DDG-5	1.009 ± 0.036
DDG-6	0.972 ± 0.072

Results

The sulphur content of dried distillers grains could be analysed with the rapid CS cube with a very good precision without any limitations. The fast analysis and easy sample preparation enables the fully automated analyses of up to 60 samples within 5 hours with the rapid S cube.